
Basic Laboratory Information System (BLIS)

Background

- To address the need for accurate and reliable laboratory data, the Centers for Disease Control and Prevention (CDC) and the Computing 4 Good @ Georgia Institute of Technology group (C4G) initiated a collaboration in 2009 to develop a robust Basic Laboratory Information System (BLIS) designed specifically for the needs of service-delivery level laboratories in resource-limited settings, which represent the majority of the nearly 2,000 laboratories that are supported by PEPFAR.
- BLIS has been developed with ongoing input from CDC field staff, implementing partners, and laboratorians in Botswana, Cameroon, Ghana, Kenya, Nigeria, Tanzania, Uganda, and Zambia.
- The system aims to address the specific needs of laboratory data collection and management from specimen receipt to results reporting.

Laboratory specimen flow





 Accurate and reliable clinical laboratory test results are a critical component of a public health approach to disease management in resource-limited settings. Laboratory data are essential for clinicians to accurately assess the status of patients' health, make accurate diagnoses, formulate treatment plans, and subsequently monitor the effects of treatment. However, in the majority of PEPFAR-supported laboratories, clinical laboratory data is generally not stored in a manner in which it can be easily accessed, summarized, or analyzed. Moreover, these data are generally recorded in multiple, nonstandardized, and often somewhat illegible log books, which presents a significant challenge to any efforts to conduct analyses or reporting.



Figure 1. Specimen labeling and workflow




Figure 2. Non-standardized laboratory data registers



Figure 3. Data storage, retrieval, and analysis

- BLIS is designed to be highly-customizable by a laboratory manager and extremely user-friendly for all laboratorians.
- It is built with a MySQL back-end with a PHP front-end, which enables it to be run from a standard internet browser.
- The system can be run locally on a single machine or on a local area network. It can also be run in a hosted or wide area network configuration.

| Basic | Labora | atory Info | rmation Sys | stem | | Logged in as: testadmin | Edit Profile Logout Comments? |
|-------------------------|----------|-------------------------|--------------------|---|-----------|-------------------------|-----------------------------------|
| Home | Lab Co | nfigurations | Test Catalog | Reports | | | |
| New Lai | b Config | uration Cano | cel | | | | Time |
| Facility Name * | | RHINO | | Tips | | | |
| Location * | | Guanajuato, Mex | ico | Please select specimen types, test types and technician | | | |
| Lab Admin * testadmin 💠 | | | | | | | accounts that are handled at |
| | | | | | | | the lab facility. |
| Specimen | Types » | C | | | | | |
| Test Type | s | Specimen Types View All | | | | | |
| Technician | ns | Enter here | | | | | |
| Custom F | ields | | /hole Blood EDTA x | Plasma EDTA | x Urine x | | |
| | | Stool x Sp | utum × | | | | |
| Submit | Cancel | | | | | I | |

Basic Laboratory Information System

Logged in as: testadmin | Edit Profile | Logout | Comments?

Home Lab Configurations Test Catalog Reports

Lab Configuration updated | « Back to Configurations

| Facility Name | RHINO |
|---------------------|--|
| Location | Guanajuato, Mexico |
| Lab Admin | testadmin |
| Specimen Types | Serum Whole Blood EDTA Plasma EDTA Urine Stool Sputum |
| Test Types | Amylase Cholesterol Platelet Count Sodium Total Albumin Uric Acid White Blood Cell Count |
| Technician Accounts | rhino_tech [RHINO Technician] Read-Write |



Target TAT Values | Site: RHINO - Guanajuato, Mexico | Cancel




Figure 4. Configuration of new laboratory in BLIS

| ome Register P | atient/Specimen | Results | Search | Repo | | |
|---------------------|-------------------|------------|-------------|------|---------------|----------------|
| ecimen Registra | tion Session No | . 20100222 | -1 Cancel | | | |
| equired field | | | | | | |
| Specimen Type * | Whole Blood EDTA | ÷ |) | | Patient ID | 123456 |
| Tests * | ✓ Platelet Coun | t E | Uric Acid | | Addl ID | - |
| | ✓ White Blood (| Cell Count | | | <u>Name</u> | RHINO Tester |
| Specimen ID * | 23456 | | | | Gender | М |
| Additional ID | | | | | Age | 40 years |
| Lab Receipt Date * | 2010 - 02 | - 18 | | | Date of Birth | 1970 (partial) |
| Collection Date * | 2010 - 02 | - 17 | | | | |
| Collection Time * | 06 🗘 : 11 🕏 | hrs | | | | |
| Comments | | | | | | |
| Report To | Doctor/Hospital | ‡ | | | | |
| Doctor/Hospital Nam | e Dr. Hernandez | | | | | |
| | | | | | | |
| | | | | | | |

Figure 5. Patient and specimen registration in BLIS

 BLIS generates a wide variety of standard and customizable reports for use in the laboratory, for reporting to patients or clinicians, and for aggregate indicator-level reporting.






Figure 5. Example turn around time report




Figure 6. Example login screen for laboratory technician. BLIS aims to empower lab technicians and laboratory managers to improve the quality of results generated in the laboratory even as the quantity of testing increases.

Next steps

- Pilot implementations of BLIS are currently being planned in Uganda, Tanzania, Cameroon, and Ghana through the guidance of the respective Ministries of Health in each country, CDC headquarters and CDC incountry staff, and with implementing partners.
- The African Field Epidemiology Network (AFENET) will lead implementations in Uganda and Tanzania
- Global Health Systems Strengthening (GHSS) will lead implementations in Cameroon and Ghana.
- We aim to have a reliable, robust and free system ready for ongoing development ready to share within one year.

For further information

We welcome your feedback.

Please feel free to contact Mark DeZalia at mdezalia@cdc.gov with questions, suggestions or comments.

We also welcome you to try BLIS!

The online version can be accessed at: lis.cc.gatech.edu

Generic administrator login User: testadmin / Password: admin123

Generic technician login User: testlab1_tech1 / Password: tech123

If you would like a custom account, please contact Mark DeZalia.